

## POWER ON Use Case (UC1)

**Primary Actor:** CES Users

**Precondition:** The device is turned off

**Main success scenario:**

1. User turns on the device power switch
2. The screen lights up and the device starts normally

**Postcondition:** The device has been turned on and is waiting for user action

## POWER OFF Use Case (UC2)

**Primary Actor:** CES Users

**Precondition:** The device is turned on

**Main success scenario:**

1. User turns off the device power switch
2. The device powers off the display.
3. Device shutdown and stop all services and output

**Postcondition:** The device has been turned off and is not consuming power

## Set Timer Use Case (UC3)

**Primary Actor:** CES Users

**Precondition:** The device has been turned on but treatment has not started

**Main success scenario:**

1. User presses the clock button; the set time to cycle between 20 minutes, 40 minutes and 60 minutes.
2. Repeat step 1 until the desired time is reached.
3. User places the device on skin. Include Activate Electrodes use case.

**Postcondition:** The treatment time of this device been adjusted to the user's desired time.  
The device displays the time set by the user

## Set Waveform Use Case (UC4)

**Primary Actor:** CES Users

**Precondition:** The device has been turned on but treatment has not started

**Main success scenario:**

1. User selected comboBox of waveform selection
2. The user selects from Alpha, Beta and Gamma waveforms with the mouse

**Postcondition:** The waveform setting of this device been adjusted to the user's desired waveform.

## Set Frequency Use Case (UC5)

**Primary Actor:** CES Users

**Precondition:** The device has been turned on but treatment has not started

**Main success scenario:**

1. User selected comboBox of frequency selection
2. The user selects from 0.5Hz, 77Hz and 100Hz with the mouse

**Postcondition:** The frequency setting of this device been adjusted to the user's desired frequency.

## Set Current Use Case (UC6)

**Primary Actor:** CES Users

**Precondition:** The treatment has been started

**Main success scenario:**

1. User presses the up arrow button; the output current will increase by 100uA each time
2. Repeat step 1 until the desired current is reached.

**Postcondition:** The output current of this device been adjusted to the user's desired current.

**Extensions:**

- 1a. User presses the down arrow button.
  - 1a1. The output current will decrease by 50uA each time. Return to step 2.

**Related Information:**

- Use case can be ended at any time if the user disconnect the skin connection
- The output current range of the device is 0-500uA

## Skin Connect Use Case (UC7)

**Primary Actor:** Device Electrode

**Precondition:** The user has set the output waveform, frequency and treatment time

**Main success scenario:**

1. The user connects the electrode to the specified location on the skin
2. Electrode informs the device to display and advance the Time of Therapy.
3. Step 1 and 2 repeat until the device's allotted therapy time has completed.
4. Electrode is no longer active/does not sense skin

**Postcondition:** The device automatically shut down after end treatment

**Extensions:**

- 1a. The electrode no longer senses skin.
  - 1a1. Electrode stops outputting an electrical current.
  - 1a2. Device pauses the Time of Therapy display.
  - 1a3. If electrode senses skin once again, return to step 1. Otherwise, 5 seconds after losing the skin connection, end use case and the treatment stop.

## View Saved Therapies Use Case (UC8)

**Primary Actor:** CES Users

**Precondition:** The device has been turned on

**Main success scenario:** User view control panel treatment history table

**Postcondition:** The user can view the saved treatment records, including the start time, treatment duration, output current, frequency and waveform

## Clear Saved Therapies Use Case (UC9)

**Primary Actor:** CES Users

**Precondition:** The device has been turned on

**Main success scenario:**

1. The user clicks the clear history button
2. The device deletes all currently saved therapies.

**Postcondition:** All currently saved therapies are deleted.

## Save a Therapy Use Case (UC10)

**Primary Actor:** CES Users

**Precondition:** Treatment is in progress

**Main success scenario:**

1. User presses the Lock button
2. The device saves the current therapy settings to the database

**Postcondition:** A desired therapy has been saved to the device's History of Treatments.

**Related Information:** The user can only record once in a treatment process, and it is forbidden to record repeatedly

## Battery undervoltage Use Case (UC11)

**Primary Actor:** Battery sensor

**Precondition:** Low battery

**Main success scenario:**

1. When device battery is less than 5 percent, system displays a 5 percent low battery icon
2. When device battery is less than 2 percent, system displays a 2 percent low battery icon and flash. After a short while, the device will be automatic shutdown.

**Postcondition:** User replenishes battery

**Related Information:**

When the battery power is higher than 5 percent, the low battery icon is not displayed

## Device Error Use Case (UC12)

**Primary Actor:** Current sensor

**Precondition:** Equipment failure causes the current to be greater than 700uA

**Main success scenario:**

1. The user simulates a device failure and presses the Error trigger button
2. The device is permanently shut down and displays a fault icon

**Postcondition:** User can no longer use this device (Need to restart the emulator to recovery)

## Not Use Timeout Use Case (UC13)

**Primary Actor:** System timer

**Precondition:** The user turns on the power and is not in the process of treatment

**Main success scenario:**

1. No operation for 30 minutes (Replace time with seconds during simulation)
2. Automatic shutdown of the device

**Postcondition:** --

**Related Information:** When the user presses any button, it will exit this use case.

